



## DIRECT MERCURY DETERMINATION IN CRUDE OIL AND OIL PRODUCTS

ASTM D7622-10(2015)  
SN/T 4429.2-2016

### INTRODUCTION

The mercury concentration in crude oil and petroleum products can vary in a wide range of less than 0.1 ppb to dozens ppm. Direct mercury determination in crude oil and petroleum products at the range above 5 ppb is covered by **ASTM D7622-10(2015) "Standard test method for total mercury in crude oil using combustion and direct cold vapor atomic absorption method with Zeeman background correction"** and **SN/T 4429.2-2016 "Determination of total mercury in crude oil – Zeeman correction – Cold atomic absorption spectrometry"**. The complex organic matrix impedes conventional quantitative analysis for mercury.

### MEASUREMENT METHOD

The Lumex Instruments method specifies a practical guidance for the standard method implementation for crude oil and petroleum products analysis using Lumex Instruments mercury analyzers.

The principle of the method is based on the reduction of Hg(II) to the atomic state due to the thermal decomposition of the mercury compounds and the follow-up transporting of mercury atoms into the analytical cell of the analyzer by the air flow. The mercury concentration is then determined from the absorption of the 254-nm resonance radiation by mercury atoms measured by the **RA-915M** or **RA-915+ mercury analyzer** combined with **PYRO-915+ attachment** or by **RA-915F mercury analyzer** using differential atomic absorption spectroscopy with Zeeman correction for background absorption.

### ANALYTICAL CHARACTERISTICS

<b>Sample composition</b>		crude oil, condensate, naphtha, gasoline, diesel fuel, lubricants, etc.
<b>Sample mass</b>		20–150 mg
<b>Measurement range</b>	<b>ASTM D7622-10 (2015)</b>	from 5 to 350 ppb
	<b>Lumex Instruments method</b>	from 5.0 to 10000 ppb
<b>Measurement time</b>		1–5 min

### ANALYSIS FEATURES

- Direct rapid analysis (1–5 min).
- Wide dynamic measurement range, no "memory effect".
- The SRM of any composition can be used for calibration and QA/QC.
- Control of the non-selective absorption during the measurement excludes analysis errors.
- No need for cylinders with compressed oxygen or other carrier gas.
- Low running cost.

### EQUIPMENT AND REAGENTS

The following equipment and materials are used for analysis:

- Mercury analyzer RA-915M or RA-915+ combined with PYRO-915+ attachment or RA-915F mercury analyzer;
- PC with Windows® XP/Vista/7/8/10 and RAPID software;
- Lumex Instruments kit, order **No 0300003044**.




## EXAMPLES OF ANALYSIS

### Direct measurements of the mercury concentration in crude oil and naphtha

Samples	C <sub>Hg</sub> , ppb	C <sub>Hg</sub> (av.), ppb	RSD, %	Samples	C <sub>Hg</sub> , ppb	C <sub>Hg</sub> (av.), ppb	RSD, %
Crude 01	760 769 752	760	0.9	Crude 06	56.9 54.8 53.9	55.2	2.3
Crude 02	19.2 16.4 17.3	17.6	6.6	Crude 07	20.4 21.8 19.9	20.7	3.9
Crude 03	8.5 7.8 7.0	7.8	7.8	Crude 08	82.3 85.0 83.9	83.7	1.3
Crude 04	74.9 68.6 65.9	69.8	5.4	Crude 09	159.5 161.9 163.2	161.5	0.9
Crude 05	1.5 1.7 1.2	1.5	14	Naphtha 115-135	21.1 21.5 22.3	21.6	2.3

### Direct mercury determination in sample Crude 01



The contents of this paper are subject to change without notice.

To get more specific information, please contact the representative at [sales@lumexinstruments.com](mailto:sales@lumexinstruments.com)